



DOCTORS
OPTICAL LABS

Featuring Advanced Polarized Technology:













- Blocks horizontal glare, boosting visual clarity.
- Blocks harmful UVA and UVB light rays.
- Backside UV-blocking anti-reflective coating reduces eye strain.
- Increases contrast with minimal color distortion.

SUN
PACKAGES




OFFERING A COMPREHENSIVE SUN PORTFOLIO



MIRROR OPTIONS:




 Black Availability: Polycarbonate Base Color: Gray 92% Blue Light Filtered	 Chrome Availability: Polycarbonate Base Color: Gray 89% Blue Light Filtered	 Blue Availability: Polycarbonate Base Color: Gray 91% Blue Light Filtered	 Gold Availability: Polycarbonate Base Color: Brown 94% Blue Light Filtered
 Purple Availability: Polycarbonate Base Color: Gray 96% Blue Light Filtered	 Red Availability: Polycarbonate Base Color: Brown 96% Blue Light Filtered	 Green Availability: Polycarbonate Base Color: Gray 89% Blue Light Filtered	 Rose Gold Availability: Polycarbonate Base Color: Brown 97% Blue Light Filtered
 Ice Blue Availability: Polycarbonate Base Color: Gray 95% Blue Light Filtered	 Pink Availability: Polycarbonate Base Color: Gray 90% Blue Light Filtered	 Silver Availability: Polycarbonate Base Color: Gray 90% Blue Light Filtered	 Champagne Availability: Polycarbonate Base Color: Gray 87% Blue Light Filtered

MATTE MIRRORS:



 Sunny Mist Availability: Polycarbonate Base Color: Brown 97% Blue Light Filtered
 Sea Mist Availability: Polycarbonate Base Color: Brown 88% Blue Light Filtered
 Charcoal Mist Availability: Polycarbonate Base Color: Gray 88% Blue Light Filtered

Pre-mirrored lenses allow for faster turn around time

SOLID COLOR OPTIONS:

 Gray 3 Availability: Polycarbonate 86% Blue Light Filtered	 Brown 3 Availability: Polycarbonate 96% Blue Light Filtered	 Gray-Green Availability: Polycarbonate 89% Blue Light Filtered
--	---	---

GRADIENT OPTIONS:

 Gray Gradient Availability: Polycarbonate 82% Blue Light Filtered	 Brown Gradient Availability: Polycarbonate 95% Blue Light Filtered
---	--

* Gold, Red, Rose Gold, Sunny Mist and Sea Mist are only available on a brown polar lens.

Pricing is for polycarbonate with backside anti-reflective coating.	Mirror & Matte Mirror		Gradient Colors (No Mirror)		Solid Colors (No Mirror)	
	List	Max Discount	List	Max Discount	List	Max Discount
EyeFOCAL OMNI NAL®	\$303.00	\$242.40	\$303.00	\$242.40	\$258.00	\$206.40
EyeFOCAL FX	\$293.00	\$234.40	\$293.00	\$234.40	\$248.00	\$198.40
EyeFOCAL SBA	\$282.00	\$225.60	\$282.00	\$225.60	\$230.00	\$184.00
EyeFOCAL SELECT+	\$271.00	\$216.80	\$271.00	\$216.80	\$232.00	\$185.60
EyeFOCAL ANTI-FATIGUE	\$241.00	\$192.80	\$241.00	\$192.80	\$184.00	\$147.20
EyeFOCAL SINGLE VISION	\$214.00	\$171.20	\$214.00	\$171.20	\$158.00	\$126.40


Requirements:

- Minimum fitting height 14mm with minimum of 10mm above fitting cross.
- All EyeFOCAL lenses come with compensated Rx.
- Pricing is for Polycarbonate and includes backside anti-reflective coating.
- For additional charges like wrap frames, refer to the National Price List.

Rx range:

- -12.00D to +6.50D, max cyl -6.00D
- Adds +0.75D to +3.50D (Anti-Fatigue +0.50D to +1.25D)

NAL is a USPTO registered trademark by Quest Vision Care Specialty Lab.



702.463.6373
 doctorsopticallabs.com
 7330 Eastgate Road, Suite #130, Henderson, NV 89011